



Disobliteration System

Reperfusion

Improving Outcomes

Better Economics

A percutaneous, one-size fits all thrombectomy/embolectomy system designed to replace surgical and drug therapies

Ideal for acute limb ischemia, organized thrombus, thrombosed grafts, embolic events, and post atherectomy debris removal



Features & Benefits

7F System

- Non surgical approach eliminates general anesthesia, lengthy ICU time, and surgical complications
- Percutaneous approach performed under Fluoroscopy provides clear endpoints
- System promotes compatibility between sheath, guidecatheter, and rapid exchange balloon

One size fits all (peripheral)

- System designed to treat above the knee and below
- Rx balloon can treat up to 10mm vessels

Control Sheath

- 7F Sheath with 10mm occlusion balloon controls flow which promotes a reduction in distal embolization
- Maintains access to treat potential underlying lesions

Disobliteration (disō-blitē-rāshōn)

Removal of a blood clot or atherosclerotic plaque from an obstructed blood vessel.

Atraumatic mechanism of action

- Balloon sweep is a proven means of clot removal and proven gentle on vessels
- Eliminates large volume blood loss associated with aspiration
- Reduces or eliminates drug therapy which potentiates a reduction of common bleeding complications associated with lytics
- Further promotes a reduction in distal embolization common with aspiration and ultrasonic therapies
- Avoids tubular necrosis, acute renal insufficiency, and potential renal failure often associated with aspiration, maceration technologies, or anything moving blood too rapidly

No need for capital equipment

- Eliminates added cost, maintenance, and eliminates connecting tubing and other accessories

Elegant catheter based design

- Easily adoptable by interventionalists
- Rx balloon nitinol core allows for tremendous pushability

Order today

314.787.7475


www.euphratesvascular.com



Product Details

Included with ICHOR Reperfusion System (ICH7FR)

1. 7F x45cm Control Sheath with 10mm Occlusion Balloon
2. 7Fx90cm Guide Catheter with 10mm Expandable Nitinol Funnel
3. 014 Rapid Exchange Catheter with 10mm Balloon and Nitinol Core
 - 035 Sheath Dilator
 - (2) Hemostatic Valves w/ side port
 - (2) 1cc Syringes



Ordering Information

ICHOR Reperfusion System..... ICH7FR
 Quantity5/bx

Customer Service: 314.787.7475
 20 South Sarah Street
 St Louis, MO 63108
www.euphratesvascular.com